



World
Leading
Manufacturer



ACTIVESHape


VACCUM ROLL MASSAGE MACHINE

WWW.VACUACTIVUS.COM

What is vacuum roller massage?

Vacuum roller massage is a cellulite reduction treatment, which was first developed over 30 years ago. Vacuum roller massage Treatment (also known as Electronic Lymphatic Drainage): Anti-cellulite, reduce fluid retention, improve skin texture, increase lymphatic and blood circulation, reduce muscle tension. Roller massage concept was found to help reduce the appearance of scar tissue. But after establishing a many health care experts realized that Vacuum roller massage technology is very effective for cellulite reduction and improve skin elasticity. Vacuum roller massage treatment involves using motorized rollers to massage and manipulate the skin tissue. A benefit of Vacuum roller massage is a non-surgical treatment as called as an one of non-invasive liposuction method. Anti-cellulite treatment is not the only benefit that can be reached using Vacuum roller massage machine ActiveShape. ActiveShape session can also help alleviate minor muscles, improve skin condition and improve blood and lymphatic circulation on treated area.

This modern technique allows to reduce the appearance of cellulite by activating natural processes of fat tissue drainage. Using head equipped with independent and motorized rolls, the skin finds its smooth and flexible aspect. ActiveShape vacuum roller massage is a non-invasive treatment to treat localized fat pockets and shape the body.



How this cellulite reduction technology works?

Cellulite is the expression used to describe pockets of fat which are trapped and cause dimpling in the skin. This dimpling is irregular and patchy and has been identified with orange peel. 90% of post-adolescent women go on to develop cellulite at some time during their life. The thighs, buttocks and the abdomen are the most common areas for a woman to develop cellulite. Although it is often mistaken for obesity, cellulite is not actually obesity related because it can also take place in thin lean women. Cellulite tends to develop in the subcutaneous fat layers. This layer of fat is unique in its structure compared to the other layers because its fatty parts are structured into specific chambers by strands of linked tissue around it. Hormones also help out with the storage of fat along with the body's metabolism in the subcutaneous fat layer. Cellulite is for the most part found in women during their hormonal changes such as puberty, pregnancy, and menopause and also whilst taking birth control pills. The vacuum roller helps to eliminate the metabolic toxins caused by intense effort to help prevent delayed onset

muscle soreness . Softening muscle tissues the rollers increase fluidity and precision of movement.

Vacuum roller therapy promotes weight loss and full-body toning by stimulating targeted areas and running them along the lymphatic lines. These special areas provide your circulatory system with the energy required to metabolize body fat, including the fat in your problem areas and the fat around your internal organs. A more efficient circulation can facilitate the removal of toxins and excess water that are causing or increasing your weight problems.

The suction created helps to reach fat located in the deeper level, and separate it from the fibers.



Active Shape vacuum roller massage benefits:

- Promote the production of collagen and elastin
- Excess water retention
- Help to fade away scars, stretch marks and orange peel
- Reduce & Prevent Cellulite

- Improve Skin Texture & Elasticity
- Eliminate Toxins
- Stimulate Blood & Lymphatic Flow
- Tone & Re-contour Body Shape
- Weight lost benefits
- No injections, anesthesia, or probing
- Massages can relieve stress and tension, and improve circulation.
- Gives it a tight, toned appearance, improves blood flow
- No downtime or recovery period
- Painless
- Visible results after 8-10 sessions (results may vary)



What is Vacuum Roller Massage Body Treatment procedure?

Vacuum roller massage is a relaxing treatment performed by specially trained therapists, with a hand-held, massaging head that rolls over and gently kneads and sucks the affected area and subjacent soft tissue. This unique combination of gentle suction and deep-reaching massage is

applied to the affected areas in a manner that cannot be replicated by any other treatments. During treatments, the technician will roll the head on the body of the individual wearing a special wear. The trained employee will adapt the intensity of the treatment to amount of cellulite present and will focus on areas that need a deeper treatment. Stimulated by motorized massaging rolls, fibroblasts release collagen and elastin. This allows fat cells to detach from their current position, leaving space for a smooth, rejuvenated and visibly clearer skin.

A such kind of anti-cellulite massage gently draws up a skin fold and rolls it, comfortably manipulating tissues to break up pockets of cellulite and stimulate circulation through all the skin layers. The increased blood and lymph flows bring fresh oxygen and nutrients to all the layers of the skin and even to the tissues below the skin. This improved

circulation then taken away the excess water, fats and cell wastes freed up by the massage. Customers noticed the experience as very relaxing - much like a deep body massage. The massaging motions of the motorized rollers break up and smooth trapped pockets of fat, which cause the skin to look firmer. The rollers move skin in three different ways: The rolling motion rolls the skin inside, to target and slim areas of excess fat. The movement of the ActiveShape massagers improves blood flow and taken out toxins from the body. It also stimulates the body's natural production of collagen - leading to firmer, tighter skin.



Can the treatment be dangerous?

Since vacuum roller massage is a non-invasive treatment, the side effects associated with the procedure are relatively minor. Patients may experience light bruising and broken capillaries at the treatment site, but vacuum roller massage with properly adjusted suction level is a painless and does not leave any bruises on the skin.

ActiveShape vacuum roller massage treatment duration recommendation

Manufacturer recommends a cure of at least 10-15 treatments of 30-45 minutes each with 3 treatments every week. Once the cure is completed, it is recommended to keep maintenance sessions every 3-4 weeks to maintain the results. The first month of treatment is very important as it will send signals to awaken cells. It is therefore crucial to follow the treatment prescription. Note that this treatment is always more efficient when accompanied by a good life hygiene, healthy eating habits and physical activity.




Vacuum roller massage treatment costs

Compared to surgical procedures such as liposuction, any vacuum roller massage is quite inexpensive. The cost usually runs from \$50 to \$75 per treatment session depends from your area. Customers can expect to get the best deal when they purchase multiple treatment sessions as a package.

What can I expect After treatment

The lymph system is stimulated during this treatment. It is responsible for removing waste from the body. That's why you want to make sure that you drink a lot of water after getting this treatment. You will be able to return to your normal activities right after treatment.

After Active Shape treatment, it is not uncommon to experience sore skin and what looks like bruising. This is both an expected and wanted effect of the session. The color and pattern of the marks depend on the level of stagnation in the area. If marks do appear, they can range from a bright red to dark purple, usually lasting 3 days to a week.



How to maintain the results of Vacuum roller massage?

There are many things that you can do to maintain the results of vacuum roll massage. While diet and exercise alone will not fight cellulite, the results of this procedure will last

longer if you live a healthy lifestyle. Try to avoid wearing clothing that is very tight. Restrictive clothing can reduce blood flow, which can have an adverse effect on your skin.

How Often Should I Come For Maintenance Sessions And For How Long?

With ActiveShape vacuum roll massage treatments we recommend to take approximately 20 sessions for the first package. You might need more, or you can go right on to do maintenance. It all depends on your blood circulation, water retention, your fat cells and the fibrosis of the connective tissue.

The Vacuum Roller massage recommends that you go on maintenance only when you feel that you have achieved the maximum results. If you had 2 sessions a week, we will suggest reducing the sessions gradually. For example, going from twice a week to once a week, then once every two weeks, then continue to once every three weeks and finally, once a month. You should keep coming once a month in order to maintain these great results.

What Is ActiveShape From VACUACTIVUS



World
Leading
Manufacturer

This is a newest machine from Vacuactivus.

Active Shape vacuum roller massage is a noninvasive method of vacuum therapy combined performance of the mobile rolls, and the radio wave bipolar.

The treatment is completely non-invasive device, gives great results in the fight against cellulite, accelerates metabolism, increases trophies tissues, reduces body fat and shapes your body.

Vacuumping skin-fold expands blood vessels, increases metabolism or the rate of removal of toxins from the cell.

Lymph is drained through the movement of the head roller of the treatment. Bipolar radio wave acts directly on the skin and subcutaneous tissue.

Tissue warms to 60 degrees Celsius, by contracting the collagen fibers. The skin becomes more resilient and definitely thicker.

Technical data

Technologies - RF, vacuum, self-propelled massage rollers

Under-pressure - 690Mbar

Radio frequency - 500kHz

Operating frequency - 0-20Hz

Display - LCD touch screen

Power consumption - 750W

Power supply - 110V for US/Canada/Japan and 220V for Europe, Australia and other countries

Noise level - 65dB


Dimensions - 210x85x80cm 120kg

RF, vacuum, self-propelled massage rollers 690mbar



Distribution of programs: Effect: relaxation, cellulite, firming, wellness Area: calves, thighs and buttocks, back, belly



Treatment Programs:



Ability to determine the parameters:

- time
- pressure
- frequency
- mode
- RF intensity

The choice of the treatment head and the position of the couch is done with the touch screen display.



The kit includes:

- unit main
- massage head (RF) + 1 applicator,
- massage head (massage) + 2 applicators, - dress / suit 5 pieces,
- a set of filters

How big room do you need for install it?

From 6m² (70sf)

HOW MUCH ELECTRICITY DOES THE DEVICE USE?

ActiveShape is extremely economic in terms of electricity consumption, because a newest energy saving technologies used for create this machine. It could be plug-in to any regular outlet (120V America/Japan or 240V Europe/Asia/Australia depends from version).

HOW LONG IS THE WARRANTY?

A standard warranty terms is available on VACUACTIVUS manufacturer web page using this link <https://vacuactivus.com/warranty-for-exercise-equipment>

WHERE CAN I TRY MACHINE?

Please contact sales department VACUACTIVUS via e-mail:

vacuactivus@gmail.com for request information about closest location or soonest trade show where available to try and review ActiveShape machine and other VACUACTIVUS products.



A SAMPLE OF BUSINESS PLAN FOR ActiveShape

HOW MUCH PRICE FOR ACTIVESHape TREATMENT?

Vacuum roll massage treatment prices have numerous variables depending on location, reputation and whether the treatment is carried out in combination with other weight loss or cellulite reduction treatments for example. The 'Single session' in the table below is only for basic treatments. Roll vacuum massage treatments are a support for get rid fat tissue or fighting against cellulite. Average is around USD60-100 depends from

PRICING SUPPORT AESTHETIC TREATMENTS

ActiveShape Clinic Turnover Per Month

What did it mean to business financially? Real Life figures using ActiveShape are currently experiencing showing the fantastic payback of a ActiveShape device in spas and weight loss centers setting.

	PREVIOUS	CURRENT
Number of patients treated per week	43	61
Number of patients treated per month	172	244
Average treatment price	\$60 per session	\$85 per session
Total monthly gross fee income	\$10,320	\$20,740

DEVICE PAYBACK CALCULATION

	Per month
Increase In Gross Turnover	+\$10,420
Current expenses	-\$900
Other monthly expenses due add cryotherapy	-\$5,000
Increase In Net Turnover .	\$4,520
Price of ActiveShape	\$12,999
Payback time	Less than 4 months

TECHNICAL SUPPORT INFORMATION

Technical support is available all around US (our technicians are locally based in a major city areas as Los Angeles, Chicago, Dallas, New York, Miami) as well as all around Europe including UK. We could send our technicians to any destination all around the world for maintain and install

equipment. Please contact our technical team for more information
vacuactivus@gmail.com

SAMPLE OF CALCULATION FOR AESTHETIC CLINIC

Parameters	used for	the example calculation:				A	B
Treatment	Mins of Sessions	Expenses	Other Costs	Income	Profit	Sessions Per Week	Sessions Per Week
Cellulite reduction treatment	15	\$3,90	\$1.00	\$70	\$65.10	1	2
Skin tightening	15	\$3,90	\$6.20	\$100	\$89.90	1	1
Body Countouring	20	\$5,20	\$4.00	\$75	\$65.80	1	2
Cellulite reduction							
Skin Recovery	10	\$2,60	-	\$40	\$37.40	1	1
Roller							
Treatment	10	\$2,60		\$40	\$37.40	1	2
Vacuum massage	15	\$3,90	\$0.50	\$80	\$122	1	1
						6	9

DEVICE PAYBACK CALCULATION	A	B
Treatment per week	6	9
Income per week	\$405	\$590
Costs per week	\$37.6	\$56.80
Profit per week	\$367.40	\$533.20
Profit per month (4.5 weeks)	\$1,653.30	\$2,399.40
Profit per year (52 weeks)	\$19,104.80	\$27,726.40

DEVICE PAYBACK CALCULATION	A	B
ActiveShape cost	\$12,999	\$12,999
Polar Bear Payback period	9 1/2 months	6 1/2 months

PURCHASE OPTION OF ActiveShape \$12,999

9,5% sales tax will be charged for sales inside California if Reseller certificate is not not provides, NO SALES TAX applied for any sales outside California). Any purchase within European Union is a subject of VAT or business clients need to have valid local VAT number. No extra taxes for export to any country outside California or European Union.

LEASING AND FINACE OPTION

VACUACTIVUS is a manufacturer and could provide just purchase option from our side. But our partners could offer a leasing and finance options on our products (available for US an UK based clients). Finance could have a different terms from a12 to 36 months with possibly discounted interest rate. Any leasing or finance offer is a subject of approval from Lender and depends from financial paperwork



of Lessor. Please contact sales department VACUACTIVUS vacuactivus@gmail.com to ask more information regarding leasing/finance.

In case of any extra questions, feel free to contact Sales Department of VACUACTIVUS for more information

e-mail: vacuactivus@gmail.com

Phone: +1-310 4679723 (including WhatsApp and WeChat)

European sales department:

WhatsApp +48-884-191520 / Direct phone line +48 58 731 0931

UK/Ireland +44 20 32891959

Australia +61 73 103 3833

Ukraine/Russia/Belarus +380 67 3701200

www.vacuactivus.com